

Tubes dimensions and MAXIMUM CAPACITIES

Centrifuges for general applications

Sample Volume	Dim (mm) approx.	Minicen	Microcen 24	Biocen 22	Biocen 22 R	Bioprocen 22 R	Unicen 21	Digicen 22/22 R	Consul 22/22 R	Digitor 22/22 R	Dillicten 22 R	Magnus 22 / 22 R
Microtiter plates	128x86x15/21/45	-	-	-	-	4/2 (2)	-	6/4/2 (2)	12/8/4 (2)	12/8/4 (2)	12/8/4 (2)	12/8/4 (2)
Microtiter plates (h:80mm)	128x86x15/21/45/80	-	-	-	-	-	-	-	10/6/2/2 (2)	10/6/2/2 (2)	10/6/2/2(2)	10/6/2/2 (2)
Capillaries	ø1,5x75	-	-	24	24	24	-	24	-	-	-	-
PCR strips 0,2 ml.	ø6x21	2	-	4	4	4	-	4	-	-	-	-
0,2-0,4 ml. / 0,5-0,6 ml.	ø6x45/ø8x30	16/12	10/10	32/24	32/24	72/72	28/28	72/72	144/144	144/144	144/144	144/144
1,5-2 ml.	ø11x42	12	10	24	24	72	28	72	144	144	144	144
5 ml.	ø13x75	-	12	-	8	6	32	32	72	104	168	104
5 ml. conical/ screw cap	ø17x60/68	-	-	12/6	12/6	12/6	-	12/6	-	-	-	-
5 ml. blood sample	ø13x82	-	12	-	8	6	32	32	68	104	128	104
7/10 ml. blood sample	ø13x107	-	10	-	8	6	32	32	68	104	128	104
10 ml.	ø13x100	-	10	-	8	6	32	32	72	104	168	104
10 ml. blood sample	ø16x107	-	10	-	8	6	32	32	52	80	112	80
10 ml. (hs) (1)	ø16x80	-	-	-	-	6	32	32	48	72	104	72
15 ml.	ø16x100	-	10	-	8	6	32	32	52	80	112	80
15 ml. conical	ø17x122	-	8	-	8	6	32	32	28	52	76	52
25 ml. conical	ø29x83	-	-	-	-	6	6	6	12	20	32	20
30 ml. / 30 ml. (hs) (1)	ø25x98	-	-	-	-	6	6	8	20	24	40	24
50 ml.	ø34x100	-	-	-	-	-	6	6	8	16	24	16
50 ml. conical	ø29x117	-	-	-	-	6	6	6	12	20	32	20
50 ml. (hs) (1)	ø29x108	-	-	-	-	6	6	6	12	20	32	20
80 ml.	ø44x100	-	-	-	-	-	4	4	6	8	12	8
80 ml. (hs)/ 85 ml. (hs) (1)	ø38x112	-	-	-	-	-	4	4	6	12	12	12
100 ml. / 125 ml.	ø48x100	-	-	-	-	-	4	4	4	6	12	6
200 ml.	ø60x120	-	-	-	-	-	-	-	4	6	6	6
250 ml.	ø60x135	-	-	-	-	-	-	-	4	6	6	6
400 ml.	ø80x118	-	-	-	-	-	-	-	4	4	4	4
500 ml.	ø90x120	-	-	-	-	-	-	-	-	4	4	4
750 ml.	ø96x130	-	-	-	-	-	-	-	-	4	4	4
1.000 ml.	ø111x128	-	-	-	-	-	-	-	-	-	4	-
Blood bags	(3)	-	-	-	-	-	-	-	-	4	4	4
EZ Single Cytofunnel™	0,1-0,5 ml	-	-	-	-	-	-	-	-	-	-	-
EZ Megafunnel™	up to 6 ml	-	-	-	-	-	-	-	-	-	-	-
Cyto-Clips™	-	-	-	-	-	-	-	-	-	-	-	-
9/15 ml.	ø16x107	-	-	-	-	-	-	-	-	-	-	-
Butyrometers	ø25x212	-	-	-	-	-	-	-	-	-	-	-
12,5 ml. finger	ø16x105	-	-	-	-	-	-	-	-	-	-	-
100 ml. 6" conical	ø44-46x162-167	-	-	-	-	-	-	-	-	-	-	-
100 ml. 8" conical	ø36-38x195-203	-	-	-	-	-	-	-	-	-	-	-
100 ml. pear	ø58-59x157-160	-	-	-	-	-	-	-	-	-	-	-

(1) High speed tubes.
(2) Allows different configurations depending on the microplates height.
(3) Check the bags features.

Tubes dimensions and MAXIMUM CAPACITIES

Centrifuges for special applications

Sample Volume	Dim (mm) approx.	Series Digitor 22 C	Lacter 21	Plasma 22	Cyto 22	Vetcen
Microtiter plates	128x86x15/21/45	-	-	-	-	-
Microtiter plates (h:80mm)	128x86x15/21/45/80	-	-	-	-	-
Capillaries	ø1,5x75	-	-	-	-	12
PCR strips 0,2 ml.	ø6x21	-	-	-	-	-
0,2-0,4 ml. / 0,5-0,6 ml.	ø6x45/ø8x30	-	-	-	-	-
1,5-2 ml.	ø11x42	-	-	-	-	6
5 ml.	ø13x75	-	-	-	-	-
5 ml. conical/ screw cap	ø17x60/68	-	-	-	-	-
5 ml. blood sample	ø13x82	-	-	8	-	-
7/10 ml. blood sample	ø13x107	-	-	-	-	-
10 ml.	ø13x100	-	-	-	-	-
10 ml. blood sample	ø16x107	-	-	-	-	-
10 ml. (hs) (1)	ø16x80	-	-	-	-	-
15 ml.	ø16x100	-	-	-	-	-
15 ml. conical	ø17x122	-	-	-	-	-
25 ml. conical	ø29x83	-	-	-	-	-
30 ml. / 30 ml. (hs) (1)	ø25x98	-	-	-	-	-
50 ml.	ø34x100	-	-	-	-	-
50 ml. conical	ø29x117	-	-	-	-	-
50 ml. (hs) (1)	ø29x108	-	-	-	-	-
80 ml.	ø44x100	-	-	-	-	-
80 ml. (hs)/ 85 ml. (hs) (1)	ø38x112	-	-	-	-	-
100 ml. / 125 ml	ø48x100	-	-	-	-	-
200 ml.	ø60x120	-	-	-	-	-
250 ml.	ø60x135	-	-	-	-	-
400 ml.	ø80x118	-	-	-	-	-
500 ml.	ø90x120	-	-	-	-	-
750 ml.	ø96x130	-	-	-	-	-
1.000 ml.	ø110x135	-	-	-	-	-
Blood bags	(3)	-	-	-	-	-
EZ Single Cytofunnel™	0,1-0,5 ml	-	-	-	12	-
EZ Megafunnel™	up to 6 ml	-	-	-	12	-
Cyto-Clips™	-	-	-	-	12	-
9/15 ml.	ø16x107	-	-	8	-	-
Butyrometers	ø25x212	-	12	-	-	-
12,5 ml. finger	ø16x105	28	-	-	-	-
100 ml. 6" conical	ø44-46x162-167	8	-	-	-	-
100 ml. 8" conical	ø36-38x195-203	8	-	-	-	-
100 ml. pear	ø58-59x157-160	4	-	-	-	-